



青岛爱普国际贸易有限公司
ALL PRO CORPORATION

Add: No.592 Tonghe Road, Qingdao, 266035 China

Tel: 86-532-5568 9612 Fax: 86-532-5568 9611 E-mail: medical@allpro.com.cn www.allprocorporation.com



ENDOTRACHEAL TUBE Oral/Nasal

Description:

- Available both with cuff and without cuff.
- Made of transparent thermosensitive PVC that softens at body temperature.
- With high-volume, low-pressure cuff.
- Murphy eye to avoid complete respiratory obstruction.
- Radiopaque line throughout the tube for X-ray visualization.
- Package:20pcs/box, 200pcs/carton, 60*37*34cm

Type	Product No	Size(#)
Oral/Nasal Without Cuff	AP20C	2
	AP25C	2.5
	AP30C	3
	AP35C	3.5
	AP40C	4
	AP45C	4.5
	AP50C	5
	AP55C	5.5
	AP60C	6
	AP65C	6.5
	AP70C	7
	AP75C	7.5
	AP80C	8
	AP85C	8.5
	AP90C	9
	AP95C	9.5
AP10C	10	

Type	Product No	Size(#)
Oral/Nasal With Cuff		
	AP30C	3
	AP35C	3.5
	AP40C	4
	AP45C	4.5
	AP50C	5
	AP55C	5.5
	AP60C	6
	AP65C	6.5
	AP70C	7
	AP75C	7.5
	AP80C	8
	AP85C	8.5
	AP90C	9
	AP95C	9.5
AP10C	10	



青岛爱普国际贸易有限公司
ALL PRO CORPORATION

Add: No.592 Tonghe Road, Qingdao, 266035 China

Tel: 86-532-5568 9612 Fax: 86-532-5568 9611 E-mail: medical@allpro.com.cn www.allprocorporation.com



REINFORCED ENDOTRACHEAL TUBE Oral/Nasal

Description:

- Available both with cuff and without cuff.
- Made of transparent thermosensitive PVC that softens at body temperature.
- With high-volume, low-pressure cuff.
- Murphy eye to avoid complete respiratory obstruction.
- Radiopaque line throughout the tube for X-ray visualization.
- A stainless spring is inserted in the tube to minimize the risk of kinking or crushing.
- Straight reinforced endotracheal tube with preloaded stylet is very convenient for use.
- Package: 20pcs/box, 200pcs/carton, 60*37*34cm

Type	Reinforced with cuff-curved	Reinforced without cuff-curved	Reinforced with cuff-straight with preloaded stylet	Reinforced without cuff-straight with preloaded stylet	Size (#)
Both curved and straight reinforced endotracheal tubes are available. Preloaded stylet is optional.	AP30-1	AP30-2	AP30-3	AP30-4	3
	AP35-1	AP35-2	AP35-3	AP35-4	3.5
	AP40-1	AP40-2	AP40-3	AP40-4	4
	AP45-1	AP45-2	AP45-3	AP45-4	4.5
	AP50-1	AP50-2	AP50-3	AP50-4	5
	AP55-1	AP55-2	AP55-3	AP55-4	5.5
	AP60-1	AP60-2	AP60-3	AP60-4	6
	AP65-1	AP65-2	AP65-3	AP65-4	6.5
	AP70-1	AP70-2	AP70-3	AP70-4	7
	AP75-1	AP75-2	AP75-3	AP75-4	7.5
	AP80-1	AP80-2	AP80-3	AP80-4	8
	AP85-1	AP85-2	AP85-3	AP85-4	8.5
	AP90-1	AP90-2	AP90-3	AP90-4	9
	AP95-1	AP95-2	AP95-3	AP95-4	9.5
	AP10-1	AP10-2	AP10-3	AP10-4	10



青岛爱普国际贸易有限公司
ALL PRO CORPORATION

Add: No.592 Tonghe Road, Qingdao, 266035 China

Tel: 86-532-5568 9612 Fax: 86-532-5568 9611 E-mail: medical@allpro.com.cn www.allprocorporation.com



NASAL PREFORMED TUBE

Description:

Available with cuff and without cuff, and reinforced tube available

Made of medical grade transparent thermosensitive PVC.
Cuffed tube with high volume low pressure cuff.
Special formed for nasal use.

Size: 2.0/2.5/3.0/3.5/4.0/4.5/5.0/5.5/6.0/6.5/7.0/7.5
/8.0/8.5/9.0/9.5/10.0

Package: 20pcs/ctn, 200pcs/ctn, 60*37*34cm



青岛爱普国际贸易有限公司
ALL PRO CORPORATION

Add: No.592 Tonghe Road, Qingdao, 266035 China

Tel: 86-532-5568 9612 Fax: 86-532-5568 9611 E-mail: medical@allpro.com.cn www.allprocorporation.com



ORAL PREFORMED TUBE

Description:

Available with cuff and without cuff, and reinforced tube available

Made of medical grade transparent thermosensitive PVC.
Cuffed tube with high volume low pressure cuff.
Special formed for oral use.

Size: 2.0/2.5/3.0/3.5/4.0/4.5/5.0/5.5/6.0/6.5/7.0/7.5
/8.0/8.5/9.0/9.5/10.0

Package: 20pcs/box, 200pcs/ctn, 60*37*34cm



青岛爱普国际贸易有限公司
ALL PRO CORPORATION

Add: No.592 Tonghe Road, Qingdao, 266035 China

Tel: 86-532-5568 9612 Fax: 86-532-5568 9611 E-mail: medical@allpro.com.cn www.allprocorporation.com



INJECTABLE ENDOTRACHEAL TUBE

Description:

Made of medical grade PVC and latex free.

With injection port

Size: 4.5/5.0/5.5/6.0/6.5/7.0/7.5/8.0/8.5/9.0

Package: 10pcs/box, 100pcs/ctn, 45*39.5*32cm



青岛爱普国际贸易有限公司
ALL PRO CORPORATION

Add: No.592 Tonghe Road, Qingdao, 266035 China

Tel: 86-532-5568 9612 Fax: 86-532-5568 9611 E-mail: medical@allpro.com.cn www.allprocorporation.com



ENDOBONCHIAL TUBE

Description:

Three styles of bronchial cuffs are available.

Two styles of connectors, fixed and non-fixed.

Low-pressure cuffs can help to reduce the lesion of mucosa.

Available also in a set together with connector and three pieces of suction catheters.

Size: Right 26Fr/28Fr/32Fr/35Fr/37Fr/39Fr/41Fr
Left 26Fr/28Fr/32Fr/35Fr/37Fr/39Fr/41Fr

Package: 10pcs/ctn, 53*33*30cm