



青岛爱普国际贸易有限公司
ALL PRO CORPORATION

Add: No.592 Tonghe Road, Qingdao, 266035 China

Tel: 86-532-5568 9612 Fax: 86-532-5568 9611 E-mail: medical@allpro.com.cn www.allprocorporation.com



CBI Pump

Specification

Volume	60ml	100ml	150ml	200ml	275ml		
Flow velocity	0.5ml/h	1.0ml/h	2.0ml/h	3.0ml/h	4.0ml/h	5.0ml/h	6.0ml/h



青岛爱普国际贸易有限公司
ALL PRO CORPORATION

Add: No.592 Tonghe Road, Qingdao, 266035 China

Tel: 86-532-5568 9612 Fax: 86-532-5568 9611 E-mail: medical@allpro.com.cn www.allprocorporation.com



CBI+PCA Pump

Specification

Volume	60ml	100ml	150ml	200ml	275ml		
Flow velocity	0.5ml/h	1.0ml/h	2.0ml/h	3.0ml/h	4.0ml/h	5.0ml/h	6.0ml/h
DIYhandle	0.5ml/min 0.5ml/min						



青岛爱普国际贸易有限公司
ALL PRO CORPORATION

Add: No.592 Tonghe Road, Qingdao, 266035 China

Tel: 86-532-5568 9612 Fax: 86-532-5568 9611 E-mail: medical@allpro.com.cn www.allprocorporation.com



Continuous Multi-velocity CBI

Specification

Volume	60ml	100ml	150ml	200ml	275ml
Flow velocity	2, 3, 4, 5ml/h	4, 6, 8, 10ml/h	6, 8, 10, 12ml/h		