



青岛爱普国际贸易有限公司  
ALL PRO CORPORATION

Add: No.592 Tonghe Road, Qingdao, 266035 China

Tel: 86-532-5568 9612 Fax: 86-532-5568 9611 E-mail: medical@allpro.com.cn www.allprocorporation.com

ISO  CE  FDA 



## BREATHING CIRCUIT - CORRUGATED

### Description:

- Breathing Circuit, corrugated tube.
- Standard connectors (15mm, 22mm).
- Kinking resistant, very high quality.
- Size(#): Adult/Pediatric
- Remark: Length can be customized, usually 1.5m, 1.6m, 1.8m, 2.4m, 3.0m
- Package: 500pcs/carton, 65\*46\*35cm



青岛爱普国际贸易有限公司  
ALL PRO CORPORATION

Add: No.592 Tonghe Road, Qingdao, 266035 China

Tel: 86-532-5568 9612 Fax: 86-532-5568 9611 E-mail: medical@allpro.com.cn www.allprocorporation.com



## BREATHING CIRCUIT - COAXIAL

### Description:


- Coaxial breathing circuit, standard connectors (15mm, 22mm).
- Inside gas sampling line. Very high quality.
- The tube color can be customized.
  
- Size(#): Adult/Pediatric
- Remark: Length can be customized, usually 1.5m, 1.6m, 1.8m, 2.4m, 3.0m
  
- Package: 500pcs/carton, 65\*46\*35cm



青岛爱普国际贸易有限公司  
ALL PRO CORPORATION

Add: No.592 Tonghe Road, Qingdao, 266035 China

Tel: 86-532-5568 9612 Fax: 86-532-5568 9611 E-mail: medical@allpro.com.cn www.allprocorporation.com



## BREATHING CIRCUIT - SMOOTHBORE

### Description:

- Breathing Circuit, smoothbore tube.
- Standard connectors (15mm, 22mm).
- Kinking resistant, very high quality.
  
- Size(#): Adult/Pediatric
- Remark:Length can be customized, usually 1.5m, 1.6m, 1.8m, 2.4m, 3.0m
  
- Package:500pcs/carton, 65\*46\*35cm



青岛爱普国际贸易有限公司  
ALL PRO CORPORATION

Add: No.592 Tonghe Road, Qingdao, 266035 China

Tel: 86-532-5568 9612 Fax: 86-532-5568 9611 E-mail: medical@allpro.com.cn www.allprocorporation.com

ISO  CE  FDA 



## BREATHING CIRCUIT - EXPANDABLE

### Description:

- Breathing Circuit, expandable tube.
- Standard connectors (15mm, 22mm).
- Kinking resistant, very high quality.
  
- Size(#): Adult/Pediatric
- Remark: Length can be customized, usually 1.5m, 1.6m, 1.8m, 2.4m, 3.0m
  
- Package: 500pcs/carton, 65\*46\*35cm



青岛爱普国际贸易有限公司  
ALL PRO CORPORATION

Add: No.592 Tonghe Road, Qingdao, 266035 China

Tel: 86-532-5568 9612 Fax: 86-532-5568 9611 E-mail: medical@allpro.com.cn www.allprocorporation.com



## BREATHING CIRCUIT - CUSTOMIZED

### Description:

- Our breathing circuit (ventilator circuit, anesthesia circuit) can be customized at many lengths, and be equipped with Water Trap, Breathing Bag (latex or latex free), or Filter, HMEF, Catheter Mount, etc. If you have any specific requirements, please just contact us.
- We have Corrugated, Expandable, Smoothbore, Coaxial and other tubes, various types.