

Percutaneous Nephrostomy Set



Features

- ★ Common instrument in urology for sequential percutaneous expansion
- ★ Used for progressive dilatation of the nephrostomy tract
- ★ Introducer sheath allows the safety of wire placement
- ★ Smooth surface for reducing patient's trauma

Contents

- ★ Fixation Device
- ★ Dilators
- ★ Tracar Needles
- ★ Guide wires

Product No.	With Fascial Dilators	Fr
PNS-10	6F, 8F, 10F, 12F	10F
PNS-12	6F, 8F, 10F, 12F, 14F	12F
PNS-14	6F, 8F, 10F, 12F, 14F, 16F	14F
PNS-16	8F, 10F, 12F, 14F, 16F, 18F	16F
PNS-18	10F, 12F, 14F, 16F, 18F, 20F	18F
PNS-20	12F, 14F, 16F, 18F, 20F, 22F	20F