




青岛爱普国际贸易有限公司 ALL PRO CORPORATION

Add: No.592 Tonghe Road, Qingdao, 266035 China

Tel:86-532-5568 9697 Fax:86-532-5568 9611 E-mail:interventional@allpro.com.cn www.allprocorporation.com



PTCA Balloon Dilatation Catheter CE 0197



- 0.017" small tip profile for excellent crossability.
- Soft tapered tip, provides outstanding ability to encircled the guidewire, minimizing vessel injury.
Stable optimized lubricious hydrophilic coating for smooth delivery.
- Transparent hub is convenient to observe the air bubble, reducing the operation risk.
- The excellent smooth transition between hub and proximal shaft is designed for enhanced pushability and torque transfer.

Technical Parameters

| | | | |
|--|-----------------|-----------|-------|
| Catheter Design | Rapid Exchange | | |
| Balloon Compliance | Semi-compliant | | |
| Marker | 2 | | |
| Material | Nylon 12 | | |
| Nominal Pressure | 1.5~3.0mm 8atm | 3.5~4.0mm | 6atm |
| Rated Burst Pressure | 1.5~3.5mm 16atm | 4.0mm | 14atm |
| Distal Shaft (outer diameter) | 1.5~2.0mm 2.6F | 2.5~4.0mm | 2.8F |
| Proximal Shaft (outer diameter) | 1.5~2.0mm 2.1F | 2.5~4.0mm | 2.3F |
| Crossing Profile | 0.026~0.030" | | |
| Usable Catheter Length | 135cm | | |
| Minimum Guiding Catheter Compatibility | 5Fr/0.056" | | |
| Guidewire Compatible | 0.014" | | |

Ordering Information

| Length (mm) | Diameter(mm) | | | | | | | |
|----------------|--------------|----------|----------|----------|----------|----------|----------|----------|
| | 1.5 | 1.75 | 2.0 | 2.5 | 2.75 | 3.0 | 3.5 | 4.0 |
| 9 | PRX15009 | PRX17009 | PRX20009 | PRX25009 | PRX27009 | PRX30009 | PRX35009 | PRX40009 |
| 11 | PRX15011 | PRX17011 | PRX20011 | PRX25011 | PRX27011 | PRX30011 | PRX35011 | PRX40011 |
| 12 | PRX15012 | PRX17012 | PRX20012 | PRX25012 | PRX27012 | PRX30012 | PRX35012 | PRX40012 |
| 14 | PRX15014 | PRX17014 | PRX20014 | PRX25014 | PRX27014 | PRX30014 | PRX35014 | PRX40014 |
| 16 | PRX15016 | PRX17016 | PRX20016 | PRX25016 | PRX27016 | PRX30016 | PRX35016 | PRX40016 |
| 17 | PRX15017 | PRX17017 | PRX20017 | PRX25017 | PRX27017 | PRX30017 | PRX35017 | PRX40017 |
| 18 | PRX15018 | PRX17018 | PRX20018 | PRX25018 | PRX27018 | PRX30018 | PRX35018 | PRX40018 |
| 20 | PRX15020 | PRX17020 | PRX20020 | PRX25020 | PRX27020 | PRX30020 | PRX35020 | PRX40020 |
| 22 | PRX15022 | PRX17022 | PRX20022 | PRX25022 | PRX27022 | PRX30022 | PRX35022 | PRX40022 |
| 23 | PRX15023 | PRX17023 | PRX20023 | PRX25023 | PRX27023 | PRX30023 | PRX35023 | PRX40023 |
| 24 | PRX15024 | PRX17024 | PRX20024 | PRX25024 | PRX27024 | PRX30024 | PRX35024 | PRX40024 |
| 26 | PRX15026 | PRX17026 | PRX20026 | PRX25026 | PRX27026 | PRX30026 | PRX35026 | PRX40026 |
| 30 | PRX15030 | PRX17030 | PRX20030 | PRX25030 | PRX27030 | PRX30030 | PRX35030 | PRX40030 |
| 31 | PRX15031 | PRX17031 | PRX20031 | PRX25031 | PRX27031 | PRX30031 | PRX35031 | PRX40031 |
| 32 | PRX15032 | PRX17032 | PRX20032 | PRX25032 | PRX27032 | PRX30032 | PRX35032 | PRX40032 |
| 34 | PRX15034 | PRX17034 | PRX20034 | PRX25034 | PRX27034 | PRX30034 | PRX35034 | PRX40034 |
| 38 | PRX15038 | PRX17038 | PRX20038 | PRX25038 | PRX27038 | PRX30038 | PRX35038 | PRX40038 |