



青岛爱普国际贸易有限公司  
ALL PRO CORPORATION

Add: No.592 Tonghe Road, Qingdao, 266035 China

Tel: 86-532-5568 9612 Fax: 86-532-5568 9611 E-mail: medical@allpro.com.cn www.allprocorporation.com



## NASAL OXYGEN CANNULA

### Description:

- Nasal Cannula, made of clear, medical-grade PVC.
- Diameter 6mm, kinking resistant, high quality.
- Size(#):

Adult straight	Adult soft curved tip
Child straight	Child soft curved tip
Infant straight	Infant soft curved tip
Neonate straight	Neonate soft curved tip
- Package: 100pcs/carton, 35\*37\*34cm