



# 青岛爱普国际贸易有限公司 ALL PRO CORPORATION

Add: No.592 Tonghe Road, Qingdao, 266035 China

Tel: 86-532-5568 9614 Fax: 86-532-5568 9611 E-mail: laboratory@allpro.com.cn www.allprocorporation.com



www.allprocorporation.com

## AP-733+ Semi-auto Biochemistry Analyzer



### FUNCTIONS

- ◆ Stream line design, nice appearance
- ◆ Self calibration and trouble alarm function
- ◆ Absolutely open system: both of programming and reagent
- ◆ Warning tone for finished sipping and any other ready function
- ◆ Function sketches on main menu, the software is simplicity, easy and fast operation.
- ◆ Real time curve monitoring enhances the reliability of result.
- ◆ Memory for 100 programs and 20000 test results, results tracking by date.
- ◆ On-board Quality control, control results available either memory or display.

### SPECIFICATIONS

- ◆ Measuring Mode: Kinetics, End point, Fixed time, Differential, Absorbance,
- ◆ Multi-standard
- ◆ Light source: 6v/10w Halogen Lamp, life Span>3000 hours
- ◆ Standard filter: 340nm, 405nm, 505nm, 546nm, 578nm, 620nm, 670nm and one free position.
- ◆ Band width:  $\leq 8\text{nm}$
- ◆ Linearity range: 0.000-3.000Abs
- ◆ Aspiration volume: 200-1400ul
- ◆ Flow cell: 32ul, Optical diameter: 1cm
- ◆ Cross contamination rate:  $\leq 1\%$
- ◆ Temperature control: 25°C, 30°C, 37°C,  $\pm 0.1^\circ\text{C}$  and be controlled by peliter
- ◆ Language: English or any language on request
- ◆ Output: Built-in printer, RS232 for additional PC
- ◆ Display: 240\*128LCD
- ◆ Power supply: AC220V $\pm 10\%$ , 50/60Hz
- ◆ Dimension: 430(L)mm $\times$ 320(W)mm $\times$ 220(H)mm
- ◆ Net weight: 9.0kg

ALL PRO CORPORATION

Add: Rm212, No.592 Tonghe Road, Qingdao, 266035 China

Tel: (86-532) 8281 2277 Fax: (86-532) 8281 2299 E-mail: laboratory@allpro.com.cn