



ALL PRO CORPORATION

Add: No.592 Tonghe Road, Qingdao, 266035 China

Tel:86-186 5321 0180 Fax:86-532-5568 9690 E-mail:dancy@allprocorporation.com www.allprocorporation.com



K1B-1L household healthcare

Performance

Power : ≤100 W
Voltage: ~220V±22V/50Hz±1Hz
Size: 375×283×283 (L×W×H) mm
Weight: 6.0 kg
Outlet: 1
Output Pressure: 0.04-0.06 MPa
Oxygen Flow Rate: 1L/2L/3L/ 5L adjustable
Oxygen concentration: 45%-90%
Flow Accuracy:±10%
Noise: ≤45dB (A)

Functions:

Basic types:

1. Running time setting
2. Far infrared remote control
3. Oxygen inhalation
4. density negative ions functions

Nebulization Type:

Basic type functions + medicine solution nebulization by output oxygen.



AK5W Medical Oxygen Concentrator

Performance

Power : 280W in average, less than 350W in highest level
Size : 305×300×535(L×W×H)mm
Weight: 14.5 kg
Outlet : 1
Output Pressure : 0.04-0.06 MPa
Oxygen Flow Rate : (0~5)L/min. Within (0~5)L/min,
the oxygen concentration ≥90%
Noise:≤48dB (A)

Functions:

Basic Type:

1. Running time calculation.
2. Bacteria filter inside can keep away 99.999% bacteria in the air.
3. Alarms for high/low pressure, power failure, and air blockage.
4. Protection for over-heat and over-load.

Nebulization Type: Basic type functions + medicine solution nebulization by output oxygen.

Purity monitoring Type: Basic type functions + oxygen purity indication by two indicating lights on the LCD label