



青岛爱普国际贸易有限公司 ALL PRO CORPORATION

Add: No.592 Tonghe Road, Qingdao, 266035 China

Tel: 86-532-5568 9614 Fax: 86-532-5568 9611 E-mail: laboratory@allpro.com.cn www.allprocorporation.com



编号Code	描述Description	型号Type	材质Material
HSF01	采样杯 (尿) Specimen Container	30ml	杯子PS, 盖子PP/cup in PS,cap in PP
HSF02	采样杯 (大便) Stool Container	30ml	杯子PS, 盖子PP/cup in PS,cap in PP
HSF03	直口、螺口、带刻度采样杯 Screw Cap Specimen Container with Model Graduation	40ml	杯子PS, 盖子PP/cup in PS,cap in PP
HSF04	直口、螺口、无刻度采样杯 Screw Cap Specimen Container without Graduation	60ml	PP
HSF05	广口、螺口、双刻度采样杯 Screw Cap Specimen Container with Two Graduations	100ml	PP
HSF06	广口、摺盖、刻度采样杯 Press Cap Specimen Container with Graduation	40ml	PP
HSF07	广口、螺口、刻度采样杯 Screw Cap Specimen Container with Graduation	60ml	PP
HSF08	广口、螺口、刻度采样杯 Screw Cap Specimen Container with Graduation	120ml	PP
HSF11	日立杯 Sample Cups for HITACHI 7020、7150	16 × 38mm	PS
HSF12	700样品杯 Sample Cups 700,for US Beckman/Dubang	14 × 25mm,1.5ml	PS
HSF13	贝克曼样品杯 Sample Cups for US Beckman CX series	13 × 24mm,0.5ml	PS