



青岛爱普国际贸易有限公司 ALL PRO CORPORATION

Add: No.592 Tonghe Road, Qingdao, 266035 China

Tel: 86-532-5568 9614 Fax: 86-532-5568 9611 E-mail: laboratory@allpro.com.cn www.allprocorporation.com



+4°C Blood Bank Refrigerator XC-240/280L

The Blood Bank Refrigerator is a specialized refrigerating equipment for blood (whole blood) storage application, which has been used in medical industries as well as cold storage of pharmaceuticals, biological products, vaccine, and so forth. It is suitable for hospitals, blood banks and sanitation and anti-epidemic stations.

Control System:

- Microprocessor-based temperature controller, temperature range $4\pm 1^{\circ}\text{C}$, Large screen LCD display of temperature, display accuracy of 0.1°C . Room temperature ranging from 0°C to 32°C
- Automatic temperature control, Defrost automatic.
- Digital display of upper and lower temperature, two sensors.
- Perfect Audible/visual alarm system: High/low temperature alarm, Door open alarm, System failure alarm
- Power supply: 220V /50Hz 1 phase, can be changed as 220V 60HZ or 110V 50/60HZ

Structure Design:

- Upright type, Exterior made from structural plate, interior made from ABS material. Four units Caster are mounted under the bottom
- Two-layer transparent door with inert gas inside, and lockable
- Interior fluorescent lighting
- Adjustable 5-6 units shelves made of quality steel wire.
- Temperature printer is standard, Optional chart recorder, blood storage basket.

Refrigeration System:

- International Famous compressor and EBM fan motor
- Forced air circulation system
- The high-effective air-cooled condenser and the fin-type evaporator
- Refrigerant as R134a, CFC free
- Certificate: ISO9001, ISO14001, ISO13485

Specs list:

Model	Capacity	External size. (WxDxH) mm	Inside size (WxDxH) mm	Input power	Weight (N.t / G.t)
XC-240L	240 L	615*605*1615	445*525*1070	500	96/105 Kg
XC-280L	280 L	615*605*1780	445*530*1265	500	105/110 Kg

