



青岛爱普国际贸易有限公司  
ALL PRO CORPORATION

Add: No.592 Tonghe Road, Qingdao, 266035 China

Tel: 86-532-5568 9610

Fax: 86-532-5568 9611

E-mail: med@allpro.com.cn

ISO  CE  FDA 

www.allprocorporation.com



---

## ABSORBENT GAUZE BANDAGE

### Description:

1. Material: 100% cotton
  2. Width: 2.5cm, 5cm, 7.5cm, 10cm, 15cm, 20cm
  3. Length: 4.5m, 5m, 10m
  4. Density: 40's\*40's; 19\*10mesh, 19\*15mesh, 26\*18mesh, etc  
Other dimensions, width, length, and package can be produced as per requirement
  5. Description: A gauze bandage is a thin, woven fabric material that is placed over a wound to keep it clean while allowing air to penetrate and promote healing. It can be used to secure a dressing in place, or it can be used directly on a wound. These bandages are the most common type and are available in many sizes. They are good for hard to reach areas that may be difficult to dress.
-