



青岛爱普国际贸易有限公司  
ALL PRO CORPORATION

Add: No.592 Tonghe Road, Qingdao, 266035 China

Tel: 86-532-5568 9610 Fax: 86-532-5568 9611 E-mail: med@allpro.com.cn



www.allprocorporation.com



## Description:

3\*6cm pad with 5\*5cm pouch, 70% isopropyl alcohol, individually sealed swabs, 100pcs/box, 200 pcs/box, 10000pcs/carton  
Alcohol swab, with 70% IPA.  
Original solution from Exxon Mobil 99.8% solution.  
60x30mm, single pack.

### Standard packing:

English neutral packing, 100pc/box or 200pc/box,  
100box/carton;  
or as client requirement.

### Specs:

Pouch size: 50\*50±2mm  
Pouch Material: Aluminum foil paper, 73g/m<sup>2</sup>  
Non-woven size: 60\*30±2mm,  
Solution: 70% isopropyl alcohol  
Solution Weight: 0.4g- 0.5g  
The material of the box: cardboard with white face  
and mottled back; 300g/m<sup>2</sup>;