



青岛爱普国际贸易有限公司  
ALL PRO CORPORATION

Add: No.592 Tonghe Road, Qingdao, 266035 China

Tel: 86-532-5568 9610 Fax: 86-532-5568 9611 E-mail: med@allpro.com.cn



www.allprocorporation.com



## ABSORBENT GAUZE SWABS

### Description:

1. Material: 100% cotton
2. Size: 5\*5cm, 7.5\*7.5cm, 10\*10cm, 15\*15cm;  
4ply, 8ply, 12ply, 16ply;  
40's\*40's; 19\*10mesh, 19\*15mesh, 26\*18mesh, etc  
folded or unfolded  
\*Other dimensions, width, length, and package can  
be produced as per requirement
3. Description: Bleached with high temperature and high  
pressure by pure oxygen, to offer the high absorbency  
and less defective. The gauze is cut as per the designed  
size and then fold into multiple size and plies in square  
form. The gauze cutting is folded into inside of gauze  
swabs, to avoid loose threads.