



青岛爱普国际贸易有限公司
ALL PRO CORPORATION

Add: No.592 Tonghe Road, Qingdao, 266035 China

Tel: 86-532-5568 9612 Fax: 86-532-5568 9611 E-mail: medical@allpro.com.cn www.allprocorporation.com

ISO  CE  FDA 



0.3CC LO-DOSE INSULIN SYRINGE

Description:

With permanently attached needle
Latex free
Sterile
Non-toxic
Non pyrogenic
Permanently attached needle
Reduced dead space for reduced Insulin waste
Smooth plunger action for greater filling accuracy
Auto ground needle tip for greater comfort

Size: 29GX1/2" / 30GX5/16"

Package: 10pcs/bag, 100pcs/box, 500pcs/carton



青岛爱普国际贸易有限公司
ALL PRO CORPORATION

Add: No.592 Tonghe Road, Qingdao, 266035 China

Tel: 86-532-5568 9612 Fax: 86-532-5568 9611 E-mail: medical@allpro.com.cn www.allprocorporation.com



1/2CC LO-DOSE INSULIN SYRINGE

Description:

With permanently attached needle
Latex free
Sterile
Non-toxic
Non pyrogenic
Permanently attached needle
Reduced dead space for reduced Insulin waste
Smooth plunger action for greater filling accuracy
Auto ground needle tip for greater comfort

Size: 28GX1/2" / 29GX1/2" / 30GX5/16"

Package: 10pcs/bag, 100pcs/box, 500pcs/carton



青岛爱普国际贸易有限公司
ALL PRO CORPORATION

Add: No.592 Tonghe Road, Qingdao, 266035 China

Tel: 86-532-5568 9612 Fax: 86-532-5568 9611 E-mail: medical@allpro.com.cn www.allprocorporation.com



1CC LO-DOSE INSULIN SYRINGE

Description:

With permanently attached needle
Latex free
Sterile
Non-toxic
Non pyrogenic
Permanently attached needle
Reduced dead space for reduced Insulin waste
Smooth plunger action for greater filling accuracy
Auto ground needle tip for greater comfort

Size: 28GX1/2" / 29GX1/2" / 30GX5/16"

Package: 10pcs/bag, 100pcs/box, 500pcs/carton



青岛爱普国际贸易有限公司
ALL PRO CORPORATION

Add: No.592 Tonghe Road, Qingdao, 266035 China

Tel: 86-532-5568 9612 Fax: 86-532-5568 9611 E-mail: medical@allpro.com.cn www.allprocorporation.com



INSULIN PEN NEEDLES

Description:

Optimum treatment with maximum simplicity
Sterile, Non-toxic, Non-pyrogenic
Double point system
Patient point with triple beveling for easy and atraumatic penetration
Cartridge point sharpened for a perfect punch with no fragmentation
Compatible with all available pen models
Easy flow of medication without any impurities

Size: 29G(0.33)X1/2"(12mm) / 31G(0.25)X1/4"(6mm)
/ 31G(0.25)X1/3"(8mm)

Package: 10pcs/bag, 100pcs/box, 500pcs/carton