



青岛爱普国际贸易有限公司
ALL PRO CORPORATION

Add: No.592 Tonghe Road, Qingdao, 266035 China

Tel: 86-532-5568 9612 Fax: 86-532-5568 9611 E-mail: medical@allpro.com.cn www.allprocorporation.com



Pipe outside diameter			
Gauge	Colour code	O.D. mm	Length mm
26G	Brown	0.45	5/8"
25G	Orange	0.5	3/4"
24G	Purple	0.55	3/4"
23G	Blue	0.6	1"
22G	Black	0.7	1"
21G	Dark Green	0.8	11/8"
20G	Yellow	0.9	11/8"
19G	Cream White	1.1	11/4"
18G	Pink	1.2	11/4"

SCALP VEIN SET

Description:

- Butterfly shaped wings for better handling and attachment to skin
- Slip or luer set end matches well with different types of infusion sets
- Siliconed cannula for smooth penetration
- Tube length: 15-30cm
- Gauge 18g - 26g available
- ETO sterilized
- Internationally accepted color-code
- Manufacturing using automated polishing and lubrication processes
- Strict quality control procedures ensure consistent performance

Package: 100pcs/1 box, 2000pcs/carton, 56x25.5x37.5cm