



青岛爱普国际贸易有限公司  
ALL PRO CORPORATION

Add: No.592 Tonghe Road, Qingdao, 266035 China

Tel: 86-532-5568 9612 Fax: 86-532-5568 9611 E-mail: medical@allpro.com.cn www.allprocorporation.com



## NEBULIZER MASK

### Description:

Made of medical PVC.  
Without nose clip.  
Medium concentration.  
Available with anti-crush 7ft.tubing.

Size: XL/L/M/S

Package: 100pcs/ctn, 50\*37\*34cm