



青岛爱普国际贸易有限公司
ALL PRO CORPORATION

Add: No.592 Tonghe Road, Qingdao, 266035 China

Tel: 86-532-5568 9612 Fax: 86-532-5568 9611 E-mail: medical@allpro.com.cn www.allprocorporation.com



RECTAL TUBE

Description:

- Soft, frosted and kink resistant PVC tubing.
- For introduction of enema solution into rectum to release
- Atraumatic, soft, rounded, closed tip with two lateral eyes for efficient drainage
- Colour coded plain connector for easy identification of size
- Length:400mm

Item No.	Item No.(DEHP FREE)	Size(Fr)	Color Code
APF-1	F260 2400	24	Violet
APF-2	F260 2600	26	Red
APF-3	F260 2800	28	Yellow
APF-4	F260 3000	30	Black
APF-5	F260 3200	32	Blue
APF-6	F260 3400	34	Green
APF-7	F260 3600	36	Light Blue