




青岛爱普国际贸易有限公司
ALL PRO CORPORATION

Add: No.592 Tonghe Road, Qingdao, 266035 China

Tel: 86-532-5568 9612 Fax: 86-532-5568 9611 E-mail: medical@allpro.com.cn www.allprocorporation.com



SUCTION CATHETER

Description:

- Made from clear soft PVC.
- Four kind of cnnector available with T-type connector, Y-type connector, Cap-cone connector, Y-type transparent connector.
- Colour coded for size identification.
- Can With Scale line.
- Length: 500mm.

Item No.	Item No.(DEHP FREE)	Size(Fr)	Color Code
APF-1	F221 05000	5	Light Gray
APF-2	F221 06000	6	Light Green
APF-3	F250 08000	8	Blue
APF-4	F250 10000	10	Black
APF-5	F250 12000	12	White
APF-6	F250 14000	14	Green
APF-7	F250 16000	16	Orange
APF-8	F250 18000	18	Red
APF-9	F250 20000	20	Yellow
APF-10	F250 22000	22	Blue