




青島愛普國際貿易有限公司  
**ALL PRO CORPORATION**

Add: No.592 Tonghe Road, Qingdao, 266035 China

Tel: 86-532-5568 9612 Fax: 86-532-5568 9611 E-mail: medical@allpro.com.cn www.allprocorporation.com



**2-way**



**3-way**

## SILICON COATED LATEX FOLEY CATHETERS

### Description:

Made of from natural latex, silicone coated  
 Rubber valve and plastic valve for different customer needs  
 Length 400mm  
 Package: 500pcs/ctn, 52.5x41x43cm

### 2-way

Product No.	Size(#)	Color Code	Balloon Volume(mL)	Remark
AP2FC-1	12	White	5-15	Without stylet
AP2FC-2	14	Green		
AP2FC-3	16	Orange		
AP2FC-4	18	Red		
AP2FC-5	20	Yellow		
AP2FC-6	22	Purple		
AP2FC-7	24	Blue		
AP2FC-8	26	Pink		
AP2FC-9	16	Orange	30	
AP2FC-10	18	Red		
AP2FC-11	20	Yellow		
AP2FC-12	22	Purple		
AP2FC-13	24	Blue		
AP2FC-14	26	Pink		

### 3-way


Product No.	Size(#)	Color Code	Balloon Volume(mL)	Remark
AP3FC-1	16	Orange	30	3-way
AP3FC-2	18	Red		
AP3FC-3	20	Yellow		
AP3FC-4	22	Purple		
AP3FC-5	24	Blue		
AP3FC-6	26	Pink		



青岛爱普国际贸易有限公司  
**ALL PRO CORPORATION**

Add: No.592 Tonghe Road, Qingdao, 266035 China

Tel: 86-532-5568 9612 Fax: 86-532-5568 9611 E-mail: medical@allpro.com.cn www.allprocorporation.com



## ALL SILICON FOLEY CATHETERS

### Description:

Made of medical-grade silicone  
 Symmetrical balloon expands equally in all directions so as to perform its function safely and efficiently  
 Length 400mm  
 Package: 500pcs/ctn, 52x43.5x39.5cm

### 2-way

Product No.	Size(#)	Color Code	Balloon Volume(mL)	Remark
AP2FC-1	12	White	5-15	Without stylet
AP2FC-2	14	Green		
AP2FC-3	16	Orange		
AP2FC-4	18	Red		
AP2FC-5	20	Yellow		
AP2FC-6	22	Purple		
AP2FC-7	24	Blue		
AP2FC-8	26	Pink		
AP2FC-9	16	Orange	30	
AP2FC-10	18	Red		
AP2FC-11	20	Yellow		
AP2FC-12	22	Purple		
AP2FC-13	24	Blue		
AP2FC-14	26	Pink		

### 3-way

Product No.	Size(#)	Color Code	Balloon Volume(mL)	Remark
AP3FC-1	16	Orange	30	3-way
AP3FC-2	18	Red		
AP3FC-3	20	Yellow		
AP3FC-4	22	Purple		
AP3FC-5	24	Blue		