



青岛爱普国际贸易有限公司
ALL PRO CORPORATION

Add: No.592 Tonghe Road, Qingdao, 266035 China

Tel: 86-532-5568 9612 Fax: 86-532-5568 9611 E-mail: medical@allpro.com.cn www.allprocorporation.com



HEPARIN CAP

Description:

1. Used together with the arterial and venous cannula.
2. Made from fine material (medical grade PVC) and easy to puncture.
3. Latex or latex free
4. It was a tight-fitting adapter.
5. The infusion of Heparin-Sodium can prevent the reflux of blood coagulation
6. Soft colour and very smooth, without any edges and corners.
7. International luer connector, excellent On bio-compatibility.
8. ISO and CE certificate
9. Sterile blister packaging
10. Package: 4000pcs/ctn, 42*26*36cm



青岛爱普国际贸易有限公司
ALL PRO CORPORATION

Add: No.592 Tonghe Road, Qingdao, 266035 China

Tel: 86-532-5568 9612 Fax: 86-532-5568 9611 E-mail: medical@allpro.com.cn www.allprocorporation.com



HEPARIN CAP WITH CAP

Description:

1. Used together with the arterial and venous cannula.
2. Made from fine material (medical grade PVC) and easy to puncture.
3. Latex or latex free
4. It was a tight-fitting adapter.
5. The infusion of Heparin-Sodium can prevent the reflux of blood coagulation
6. Soft colour and very smooth, without any edges and corners.
7. International luer connector, excellent On bio-compatibility.
8. ISO and CE certificate
9. Sterile blister packaging
10. With cap
11. Package: 4000pcs/ctn, 42*26*36cm