



青岛爱普国际贸易有限公司
ALL PRO CORPORATION

Add: No.592 Tonghe Road, Qingdao, 266035 China

Tel: 86-532-5568 9612 Fax: 86-532-5568 9611 E-mail: medical@allpro.com.cn www.allprocorporation.com



STERILE INFUSION SETS

Description:

COMMON FEATURES

- Soft and kink resistant PVC tubing
- Clear, transparent and flexible drip chamber
- Moulded latex flush ball injection membrane for intermittent medication
- Accurate flow controller
- Approximately 15 drops / ml , 20 drops / ml, 60 drops / ml

Tube Length: 150 cm , 160 cm , 180 cm , 200 cm
Tube: I.D. : 3.0 mm , O.D. : 4.1 mm

Packing: PE package or paper-plastic package

Options available:

- With or without air vented spike
- With or without needle
- With or without "Y" injection port (Latex or Latex - free)
- Injection site material with latex or polyisoprene
- Luer lock or luer slip connector
- With or without flow regulator integrated into Infusion set for precise flow control
- With or without filter

Size: 18G/19G/20G/21G/22G/23G/24G/25G/26G/27G

Package: 500pcs/ctn, 58*38*36cm



青岛爱普国际贸易有限公司
ALL PRO CORPORATION

Add: No.592 Tonghe Road, Qingdao, 266035 China

Tel: 86-532-5568 9612 Fax: 86-532-5568 9611 E-mail: medical@allpro.com.cn www.allprocorporation.com



BLOOD TRANSFUSION SETS

Description:

Designed for transfusion of blood or blood component
Moulded latex flush ball injection membrane for intermittent medication
Approximately 20 drops / ml
Clear, transparent and soft tubing manufactured from medical grade PVC
With 200 micron filter and 18 G hypodermic needle

Packing: PE package or paper-plastic package

Options available

- With or without Airvent
- With or without Needle
- With or without "Y" site (Latex or Latex - free)
- Luer lock or Luer slip

Size: 18G/19G/20G/21G/22G/23G/24G/25G/26G/27G

Package: 400pcs/ctn, 58*40*38cm