



青岛爱普国际贸易有限公司
ALL PRO CORPORATION

Add: No.592 Tonghe Road, Qingdao, 266035 China

Tel:86-532-5568 9605 Fax:86-532-5568 9611 E-mail:diana@allprocorporation.com www.allprocorporation.com



Bouffant Cap

Disposable Operating Caps

ADVANTAGES

- Light weight, soft
- Comfortable to wear

Protect the health

ITEM	NAME	SIZE	DESCRPTTON	MATERIAL	COLOR	PACKING
1(XM1001007)	Bouffant cap with peak	18"20"21"24"	head circumference 24 g-50g/sq.m. peak 80g-100g/sq.m.	100%PP	Any colors	100pcs/polybag, 100pcs/box
2(XM1001008)	Clip cap/mob cap	18"20"21"24"	with single elastic/with double elastic, 10 g-15g/sq.m.	100%PP	Any colors	100pcs/polybag, 100pcs/box
3(XM1001009)	Sister cap	18"20"21"24"	with ties, 35g-60g/sq.m.	100%spunlace	Any colors	100pcs/polybag, 100pcs/box
4(XM1001010)	Bouffant cap (nurse cap)	18"20"21"24"	elastic, 10g-30g/sq.m.	100%PP	Any colors	100pcs/polybag, 100pcs/box
5(XM1001011)	Bouffant cap with mesh	18"20"21"24"	elastic, 10g-30g/sq.m.	Polyester	Any colors	100pcs/polybag, 100pcs/box
6(XM1001012)	Shower cap	18"20"21"24"	available in different thickness	100%PE	Any colors	100pcs/polybag, 100pcs/box