



青岛爱普国际贸易有限公司
ALL PRO CORPORATION

Add: No.592 Tonghe Road, Qingdao, 266035 China

Tel:86-532-5568 9605 Fax:86-532-5568 9611 E-mail:diana@allprocorporation.com www.allprocorporation.com



Isolation Gown

BODY WEAR

Protection Gown

- Soft and light weight material provides long term comfort for long procedures
- The comfortable cuts allow excellent freedom of movement
- Average bacterial filter efficiency 99%

ADVANTAGES

Protect the health

ITEM	NAME	SIZE	DESCRPTION	MATERIAL	COLOR	PACKING	
1(XM1004001)	PROTECTIVE GOWN	M 174X128cm	Coverall with Hood and without boots, elastic at waist and wrists	100%PPPP+PE, SMS,SF (anti-static or not anti-static)	Any colors	1pcs/PEbag	
2(XM1004002)	PROTECTIVE GOWN	2XL 189X150cm	Coverall with blue/white tape and hood, without boots, lastic cuffs	100%PPPP+PE, SMS,SF (anti-static or not anti-static)	Any colors	1pcs/PEbag	
3(XM1004003)	PROTECTIVE GOWN	XL 184X140cm	Coverall with adhesive flap and Hood ,without boots, elastic at waist and wrists	100%PPPP+PE, SMS,SF (anti-static or not anti-static)	Any colors	1pcs/PEbag	
4(XM1004004)	PROTECTIVE GOWN	L 178X132cm	cover all with collar and boots	100%PPPP+PE, SMS,SF (anti-static or not anti-static)	Any colors	1pcs/PEbag	
5(XM1004005)	PROTECTIVE SUIT	Jacket	M 66*128cm	jacket with hood and zipper elastic cuffs or knitted cuffs, elastic at the bottom of the jacket 30g-70g/sq.m.	100%PPPP+PE, SMS,SF(anti-static or not anti-static)	Any colors	1suit/PEbag
			L 70*132cm				
			XL 74*136cm				
			2XL 78*140cm				
	PROTECTIVE SUIT	Trousers	M 105*74cm	trousers with elastic on waist 30g-70g/sq.m.	100%PPPP+PE, SMS,SF(anti-static or not anti-static)		
			L 110*76cm				
			XL 115*78cm				
			2XL 120*80cm				

Remark: Spedal packing way can be provided upon request The sterilized version is available. Other different size available.