



青岛爱普国际贸易有限公司  
ALL PRO CORPORATION

Add: No.592 Tonghe Road, Qingdao, 266035 China

Tel:86-532-5568 9605 Fax:86-532-5568 9611 E-mail:diana@allprocorporation.com www.allprocorporation.com



## PATIENT GOWN

### BODY WEAR

- SMS material is comfort offer maximum protection
- Soft and light weight material provides long term comfort for long procedures
- Water-repellent, oil-repellent, alcohol-repellent, blood-repellent and antistatic
- Avoids cross infection

### ADVANTAGES



### Protect the health

ITEM	NAME	SIZE	DESCRPTTON	MATERIAL	COLOR	PACKING
1(XM1006001)	HAIRCUT SHAWL	106*138cm	with ties or magic snap on the neck PP25g-50g/sq.m, PE0.02mm-0.035mm	PP/PE/PP+PE/SMS	Any colors	1pcs/PEbag
2(XM1006002)	SAUNA DRESS / KIMONO	110*131cm	ties on the waist	PP/PE/PP+PE/SMS	Any colors	1pcs/PEbag
3(XM1006003)	PATIENT PACKS	Cap	21inch	the gown with short sleeves, round collar or "V" collar PP25g-60g/sq.m	PP/PE/PP+PE/SMS	Any colors
		Patient gown	110X130cm 115X135cm			
		Underwear	51*35cm			
		Shoe covers	38*18cm 42*18cm			
4(XM1006004)	SCRUB SUIT	Jacket	S 69*53cm	short sleeves, two pockets 40g-70g/sq.m.	PP/SMS	Any colors
			M 74*58cm			
			L 79*63cm			
		Trousers	XL 84*67cm	a pocket at the back 40g-70g/sq.m.	PP/SMS	Any colors
			S 105*55cm			
			M 110*57cm			
	L 115*59cm					
	XL 120*61cm					

Remark: Special packing way can be provided upon request. The sterilized version is available. Other different size available.