



AP-1-1

Anti decubitus mattress(lattice style)

特点:

本产品为高科技交替式换气, 超低压大流量气垫床, 专供防褥疮用。通过气垫循环波动, 使患者肌肉组织免于固定积压于同一点而产生坏死, 并能促进血液循环, 改善微循环系统, 避免人体静脉及动脉受过量压力而无法正常交换氧气及养分; 使细胞不会坏死, 最终预防和治疗褥疮。

Features:

This therapy system is designed for those with superficial established pressure sore or with existing tissue damaged, or assessed to be a low risk or be at medium risk, even at very high risk of pressure sore development. It can provide excellent pressure reduction through continuous movement. It has been proven clinically beneficial and suitable for all healthcare environment.

Specifications		AP-1-1
Pump	Electrical	AC220~240V/100~120V,50HZ/60HZ
	Air output	≥4L/min
	Pressure Range	≥12kPa(90mmHg)
	Cycle Time	5~6min
Mattress	Arrangement	2.5" bubble pad
	Material	Medical grade PVC
	Dimension(L×W×H)	190×90×6.4cm
	weight	2.6kgs



AP-1-2

Anti decubitus mattress(striped style)

特点:

本产品为高科技交替式换气, 超低压大流量气垫床, 专供防褥疮用。通过气垫循环波动, 使患者肌肉组织免于固定积压于同一点而产生坏死, 并能促进血液循环, 改善微循环系统, 避免人体静脉及动脉受过量压力而无法正常交换氧气及养分; 使细胞不会坏死, 最终预防和治疗褥疮。

Features:

This therapy system is designed for those with superficial established pressure sore or with existing tissue damaged, or assessed to be a low risk or be at medium risk, even at very high risk of pressure sore development. It can provide excellent pressure reduction through continuous movement. It has been proven clinically beneficial and suitable for all healthcare environment.

Specifications		AP-1-2	AP-1-4	AP-1-5	AP-1-6	
Pump	Electrical	AC220~240V/100~120V,50HZ/60HZ				
	Air output	≥4L/min				
	Pressure Range	≥12kPa(90mmHg)				
	Cycle Time	10~12min				
Mattress	Arrangement	22×3" cells	18×4.3" cells	17×5" cells	22×3" cells and small cells	
	Material	Cover	Nylon/PVC	Nylon/PVC	Nylon/PU	Nylon/PVC
		Cells	Nylon/PVC	Nylon/PVC	Nylon/PU	Nylon/PVC
		Base	Nylon/PVC	Nylon/PVC	Nylon/PU	Nylon/PVC
	Dimension (L×W×H)	200×90×7.6cm	200×90×10.9cm	200×90×12.7cm	200×90×7.6cm	

AP-1-3

Anti decubitus mattress(Air-jet style)


特点:

防褥疮喷气式床垫表面开有很多微孔, 利用床垫喷出的气体带走人体体表水分, 降低湿度, 从而干燥皮肤, 并置换污浊空气, 抑制细菌生长, 以及将病人体压分散在床垫上的原理, 来有效的防止褥疮的产生, 是专为减轻长期卧床病人的烦恼和痛苦, 为减轻护理人员的劳动强度而制造的。

Features:

The feature of bed sore remedy air mattress(air-jet type) is the tiny holes on the surface for decreasing humidity in bedding, and replace the foul air and control the growth of microbes, it can effectively prevent from bed sore by dispersing the weight of patients' body. It is made for decreasing the suffer and bother of long-term bedridden patients, also with lesser labor for nursing people.

Specifications		AP-1-3
Pump	Electrical	AC220~240V/100~120V,50HZ/60HZ
	Air output	≥6L/min
	Pressure Range	≥14kPa(105mmHg)
Mattress	Arrangement	Air-Jet pad
	Material	Nylon/PVC
	Dimension(L×W×H)	190×85×8.9cm



AP-1-4



AP-1-5



AP-1-6



APC-2

Cervical vertebra retractor

主要适用于颈部骨折位、滑椎、椎间盘突出、过伸过屈损伤等疾病的治疗和预防, 该牵引器凭借适度的空间压力支撑头部重量, 在固定状态下拉伸脖颈, 以此缓和颈、肩部负担, 消除酸痛。

Designed for treating and preventing neck fracture dislocation of vertebra. Protrusion of the intervertebral disc and excessive stretch and draw. It can support the patient's head efficiently depending on appropriate air pressure, so that it can release the neck and shoulder burden, dispel the pain.



APC-3-1

Anti Decubitus seat(round, square)

特点:

防褥疮座垫供久病床或下肢瘫痪的患者使用, 能有效防止褥疮生成, 采用PVC植绒料制作, 内置PVC内胆, 更具有牢固性。

Features:

It used for long-term bedridden patients or lower limbs paralysis people, and can available prevent the emergence of decubitus. This seat is made of PVC plant wool and more fastness

防褥疮座垫 Anti Decubitus seat	圆座垫round seat APC-3-1	方座垫square seat APC-3-2
膨胀尺寸 Inflate size	Φ400×Φ130mm	410×400mm
最大承载重量 The max loading weight	120kgs	120kgs



APC-3-2